Cortisol

VET control kit

For veterinary use only!

Intended use:

For quality control procedures of the quantitative determination with a Cortisol VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

| | Order information | Indication | Kit size |
|---|--|--|--|
|  | Order number: VT0290 Order number: VT0291 Order number: VT2900 | Cortisol VET test kit Cortisol VET test kit Cortisol VET control kit | 16 tests 6 tests 1 x 2 ml (decision level) |

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This Cortisol VET control kit is a valuable tool for assuring quality and precision of the cortisol determination with a Cortisol VET test kit and provides a reliable marker for monitoring the accuracy of results.

Preparation and Handling

Controls are liquid and ready to use

Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date on the label. Once opened vials are stable for 5 days when stored tightly capped at 2 - 8 °C.

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only!
DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents.

This product contains a substance with animal origin. Hazardous properties cannot be excluded. No dangers expected for the user if used according to directions. Donors have been tested and are found to be clinically healthy. Even though tests on infectious diseases ended negative, this product should be handled as potentially infectious biomaterial.

Set the species to „Control“ before running the QC material.

Compare the results with the values indicated in the lot-specific value table.

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor.

Durchführung einer Cortisol VET Kontrollmessung Processing of a Cortisol VET control measurement

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Press „Measurement“ button, enter required information using the touchscreen
- 1.3 Set species to „Control“
- 1.4 Place ERS cuvette in test kit rack
- 1.5 Place ERS cap in test kit rack
- 1.6 Take out bottle containing dilution buffer from the test kit package

2. Reagent and control preparation

- 2.1 Aspirate 400 µl from the dilution buffer bottle
- 2.2 Dispense 400 µl dilution buffer INTO THE LIQUID in the ERS cuvette

Allow single test to warm up to room temperature for at least 10 minutes before measurement! Put test kit package back into refrigerator

- 2.3 Dispense some drops of control liquid into serum cup and aspirate 40 µl using the pipette.

3. Control processing

- 3.1 Dispense 40 µl control liquid INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

Wertetabelle für das Cortisol VET Control Kit Value sheet for the Cortisol VET control kit

LOT

2307-2

Kontrollwerte | Control Values

2024-02

| Target | value (min.) | value (max.) |
|------------|--------------|--------------|
| 5.1 µg/dl | 3.1 µg/dl | 7.1 µg/dl |
| 141 nmol/l | 84 nmol/l | 197 nmol/l |

Technical details subject to change without notice.

Doc-ID: e07eaa6a-baee-4139-be32-c2b2f856c22c ; latest revision: DP 3.0 ; date of print: 2023-06-27

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