SAA VET control kit

For veterinary use only!

Intended use:

For quality control procedures of the quantitative determination with the SAA VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT0170 Order number: VT0171 Order number: VT1700	SAA VET test kit SAA VET test kit SAA VET control kit	6 tests 16 tests 1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This SAA VET control kit is intended for use as a test sample of known concentration for the SAA VET assay. The user can compare observations with expected ranges as a means of assuring consistent performance of the solo SAA VET test kit. The SAA VET control is prepared from a human serum base.

Required Material

(not supplied with the kit)

- Pipettes
- Distilled water

Preparation and Handling

Control material is lyophilized and must be reconstituted before use!

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Add exactly 2 ml of distilled/deionized water.
Warning: Inaccurate reconstitution of the control and error in assay technique can cause erroneous results.
3. Close the vial carefully and dissolve the components completely by occasional gently swirling. Avoid the formation of foam.
4. Allow the contents to reconstitute over 15 minutes before.
5. Analyse the control in the same way as a patient sample.

Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date stated on the label.

Stability after reconstitution:

at 2 - 10 °C: 2 weeks

Close immediately after use. Store control tightly capped when not in use. Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only! DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents. This control utilizes human origin serum. The serum is negative HBs antigen, HIV (HIV-1, HIV-2) and HCV antibodies test. It is impossible to completely rule out the possibility that these viruses or other pathogens may be present. Therefore when using these reagents, handle them with same degree of caution as patient serum.

This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in microbiological and biomedical laboratories" (1999).

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table


If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

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Processing of a **SAA** VET control measurement

1.

CUBE-VET SOLO

1.1: A hand holds a test strip labeled 'SAA' and inserts it into the 'SAA' slot of the 'EUROlyser' device.


1.2: A hand holds a test strip with a red cap and inserts it into a well of a multi-well plate.

1.3: A hand pours red liquid from a test strip into a well of a multi-well plate.

1.4: The device screen shows a green play button and a test strip icon. Below, a 'Measurement' button is shown with a play button icon, and a hand is shown pressing a button on the device.

2.


CUBE-VET SOLO

2.1: The device screen shows 'Specie:' with radio buttons for 'Control', 'Dog', 'Cat', and 'Young cat'. A green checkmark is next to 'Control' and a red 'X' is next to 'Cat'. Below, 'Holder:' is shown with a dropdown menu containing 'Control', 'Dog', 'Cat', and 'Young Cat'. A red checkmark is next to 'Control' and a red 'X' is next to 'Young Cat'.

2.2: A hand uses a pipette to add liquid to a test strip.


3.

3.1: A hand uses a pipette to add liquid to a test strip.

3.2: A hand holds a test strip with red liquid and inserts it into a well of a multi-well plate.

CUBE-VET SOLO

3.3: A hand inserts a test strip into the device. The screen shows 'Insert Cartridge' with a play button and a red 'X'.

3.4: A hand presses a button on the device. The screen shows a play button and a red 'X'.

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 2. Control preparation**

For further details please see user manual of laboratory photometer

 - 2.1 Set species to „Control“
 - 2.2 Aspirate 5 µl control liquid out of control vial using a pipette

- 3. Control processing**
 - 3.1 Dispense 5 µl control liquid INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply ERS cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: 2359eda6-fb10-4dd9-b62f-56c8c358b6fb ; latest revision: DP 2.0 ; date of print: 2019-10-14

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