SDMA VET control kit

For veterinary use only!

Intended use:

For quality control procedures of the quantitative determination with the SDMA VET test kit on a solo or CUBE-VET analyser

 Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

| | Order information | Indication | Kit size |
|---|--|--|--|
|  | Order number: VT0300 Order number: VT0301 Order number: VT3000 | SDMA VET test kit SDMA VET test kit SDMA VET control kit | 16 tests 6 tests 1 x 2 ml (decision level) |

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This SDMA VET control kit is a valuable tool for assuring quality and precision of the SDMA determination with the SDMA VET test kit (VT0300 / VT0301) and provides a reliable marker for monitoring the accuracy of results.

Preparation and Handling

The controls are liquid and ready to use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. To avoid contamination use single use pipetting tips only.

Stability and Storage

Unopened vials are stable at 2 – 8 °C until the expiration date stated on the label. Once opened, vials are stable for 30 days, when stored tightly capped at 2 – 8 °C. Close immediately after use. Warning: DO NOT FREEZE!

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only!
DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents.

Set the species to „Control“ before running the QC material




Compare the results with the values indicated in the lot-specific value table

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

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Processing of a **SDMA** VET control measurement

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 2. Control preparation**

For further details please see user manual of laboratory photometer

 - 2.1 Set species to „Control“
 - 2.2 Put a few drops of control liquid into a serum cup by squeezing the bottle.
 - 2.3 Aspirate 40 µl control liquid out of serum cup using a pipette

- 3. Control processing**
 - 3.1 Dispense 40 µl control liquid INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply ERS cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

Wertetabelle für das SDMA VET Control Kit Value sheet for the SDMA VET control kit

LOT

2011-1

Kontrollwerte | Control Values

2021-07

| Target | value (min) | value (max) |
|-------------|-------------|-------------|
| 60.0 µg/dl | 45.0 µg/dl | 75.0 µg/dl |
| 2.97 µmol/l | 2.23 µmol/l | 3.71 µmol/l |

Technical details subject to change without notice.

Doc-ID: 62f6f78f-a1a7-40a9-ad24-0a7b856e0315 ; latest revision: DP 3.0 ; date of print: 2020-09-24

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