For veterinary use only!

For quality control procedures of the quantitative determination with the cCRP VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT 0100 Order number: VT 0101 Order number: VT 1000	cCRP VET test kit cCRP VET test kit cCRP VET control kit	16 tests/pack 6 tests/pack 1 x 0.5 ml (high level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This cCRP VET control kit is intended for use in a veterinary clinical laboratory as a control for quantitative determination of cCRP in the high-level range. It was found to be non-reactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However this material should be handled as though capable of transmitting infectious disease.

Preparation and Handling

The controls are liquid and ready-to-use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. We recommend not to pipette directly from the original bottle as well as using single-use pipetting tips to avoid contamination. Dispense about 5 drops into a serum cup and pipette out of it.

Dispose of cup after use.

For further details on how to use the control kit please refer to the description inside this package insert.

Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date stated on the label. Once opened vials are stable until the expiry date when stored tightly capped at 2 - 8 °C.

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only!
DO NOT INGEST! Observe all necessary precautions for the use of laboratory reagents. As no testing method can completely rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table


If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

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1.

CUBE-VET SOLO


1.1: A hand holds a cCRP test strip. 1.2: The test strip is inserted into the device. 1.3: The test strip is fully inserted into the device. 1.4: The device screen shows a green play button and a measurement icon. Below the screen, a 'Measurement' button is shown with a downward arrow pointing to a smaller screen with a play button and a measurement icon.

1.1 → 1.2 → 1.3 → 1.4

2.


CUBE-VET SOLO

2.1: The device screen shows 'Specie:' with 'Control' selected (indicated by a blue dot) and 'Dog' unselected. A green checkmark is visible next to 'Control', and a red 'X' is next to 'Dog'. Below, 'Holder:' is shown with a dropdown menu containing 'Control', 'Control', and 'Dog'. A red checkmark is visible next to the first 'Control' option, and a red 'X' is next to 'Dog'. 2.2: A hand uses a pipette to add liquid to the test strip.

2.1 → 2.2

3.

3.1: A hand uses a pipette to add liquid to the test strip. 3.2: The test strip is inserted into the device. 3.3: A hand inserts a cartridge into the device. 3.4: The device screen shows 'Insert Cartridge' with a green play button and a red 'X'. Below the screen, a hand is shown inserting a cartridge into the device.

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 2. Control preparation**

For further details please see user manual of laboratory photometer

 - 2.1 Set species to „Control“
 - 2.2 Aspirate 5 µl control liquid with pipette

- 3. Control processing**
 - 3.1 Dispense 5 µl control liquid INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply ERS cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: 4b52da6a-d634-4417-9557-2d1079b40cee ; latest revision: DP 2.0 ; date of print: 2018-06-13

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