For veterinary use only!

For quality control procedures of the quantitative determination with the NH3 Ammonia VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT 0160 Order number: VT 1600	NH3 Ammonia VET test kit NH3 Ammonia VET control kit	6 tests/pack 1 x 1 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This NH3 Ammonia VET control kit is intended for use as a test sample of known concentration for the NH3 Ammonia VET assay. The user can compare observations with expected ranges as a means of assuring consistent performance of the NH3 Ammonia VET test kit.

Preparation and Handling

The controls are liquid and ready-to-use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. To avoid contamination use single use pipetting tips only.

Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date stated on the label. Once opened vials are stable until the expiry date when stored tightly capped at 2 - 8 °C.

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only!
DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents. All biological material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table


If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor


Durchführung einer NH3 Ammoniak VET Kontrollmessung Processing of a NH3 Ammonia VET control measurement

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






English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Detach cap from ERS cuvette
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 2. Control preparation**

For further details please see user manual of laboratory photometer

 - 2.1 Set species to „Control“
 - 2.2 Aspirate 40 µl control liquid using a pipette

- 3. Control processing**
 - 3.1 Dispense 40 µl control liquid INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

Wertetabelle für das NH₃ Ammoniak VET Control Kit Value sheet for the NH₃ Ammonia VET control kit

LOT 2402-1

2024-10

Kontrollwerte | Control Values

Target	value (min)	value (max)
401 µmol/l	269 µmol/l	533 µmol/l

Technical details subject to change without notice.

Doc-ID: 507ca0da-4ff6-4478-b4db-2067a5a77f89 ; latest revision: DP 2.0 ; date of print: 2018-10-03

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