## SDMA VET control kit

### Français

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### English

For veterinary use only!

**Intended use:** For quality control procedures of the quantitative determination with the SDMA VET test kit on a solo or CUBE-VET analyser.

| Order information    | Indication       | Kit size              |
|----------------------|------------------|-----------------------|
| Order number: VT0300 | SDMA test kit    | 16 tests              |
| Order number: VT0301 | SDMA test kit    | 6 tests               |
| Order number: VT3000 | SDMA control kit | 1 x 2 ml (dec. level) |

⚠ **Test kit preparation:** Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

### Summary

The use of this SDMA VET control kit serves to reliably assess the quality and precision of the SDMA determination with an SDMA VET test kit (VT0300, VT0301) and therefore enables the monitoring of the accuracy of the test results.

### Preparation and Handling

The controls are liquid and ready to use.

Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. To avoid contamination use single use pipetting tips only.

### Stability and Storage

Unopened vials are stable until the expiration date stated on the label when stored at 2 - 8 °C. Once opened, vials are stable for 30 days, when stored tightly capped at 2 – 8 °C. Close immediately after use.

Warning: DO NOT FREEZE!

### Warnings and Precautions

This control kit is for veterinary diagnostic use only!

DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents.

⚠ Set the species to „*Control*“ before running the QC material.

⚠ Compare the results with the values indicated in the lot-specific value table.

⚠ If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor.

### Waste Management

Please refer to local legal requirements

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