cCRP VET control kit

Français

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English

For veterinary use only!

Intended use: For quality control procedures of the quantitative determination with the cCRP VET test kit on a solo or CUBE-VET analyser.

Order information	Indication	Kit size
Order number: VT0100	cCRP test kit	16 tests
Order number: VT0101	cCRP test kit	6 tests
Order number: VT1000	cCRP control kit	1 x 0.5 ml (high level)

⚠ **Test kit preparation:** Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The use of this cCRP VET control kit serves to reliably assess the quality and precision of the canine CRP determination with a cCRP VET test kit (VT0100, VT0101) and therefore enables the monitoring of the accuracy of the test results.

Preparation and Handling

The control is liquid and ready to use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. To avoid contamination, use single-use pipetting tips only.

Stability and Storage

Unopened vials are stable until the expiration date stated on the label when stored at 2 - 8 °C. Once opened, vials are stable for up to 30 days, when stored tightly capped at 2 – 8 °C.

Warnings and Precautions

This control kit is for veterinary diagnostic use only!

DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents.

⚠ As no testing method can completely rule out the potential risk of infection with absolute certainty, the control material should be treated just as care fully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.

⚠ Set the species to „*Control*“ before running the QC material.

⚠ Compare the results with the values indicated in the lot-specific value table.

⚠ If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor.

Waste Management

Please refer to local legal requirements

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

1. Preparation of test system
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touch screen

2. Control preparation
For further details please see user manual of laboratory photometer
 - 2.1 Set species to „Control“
 - 2.2 Dispense some drops of control liquid into serum cup and aspirate 5 µl using a pipette.

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3. Control processing
 - 3.1 Dispense 5 µl control liquid INTO THE LIQUID in ERS cuvette.
 - 3.2 Apply ERS cap firmly onto ERS cuvette.
 - 3.3 Place ERS cartridge into laboratory photometer.
 - 3.4 Start automatic control processing by pressing the start button on the solo laboratory photometer or by closing the door of the CUBE-VET laboratory photometer

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