







scil Vet abc

Valeurs de référence



	 Chien		 Chat		 Cheval		 Bovin		 Lapin		 Furet	
	de	à	de	à	de	à	de	à	de	à	de	à
Globules blancs ($10^3/\mu\text{l}$)	6	12	5	11	5	10	7	10	6.3	10	2.5	5.5
Lymphocytes ($10^3/\mu\text{l}$)	1	3.6	1	4	1.5	4			3.3	7	0.3	1.3
Monocytes ($10^3/\mu\text{l}$)	0	0.5	0	0.5	0	0.4			0	0.4	0	0.2
Granulocytes ($10^3/\mu\text{l}$)	3	10	3	12	3	8			1.6	3.7	0.4	2
Globules rouges ($10^6/\mu\text{l}$)	5.5	8.5	5	10	5.5	11	5.5	8	5.2	6.8	6.5	12
Hémoglobine (g/dl)	15	20	8	17	9	18	10	13	11.5	15.1	15	18
Hématocrite (%)	44	57	27	47	30	50	30	40	36	47	40	60
VGM (μm^3)	60	77	40	55	37	55	45	60	65	76	45	65
TCMH (pg)	17	23	13	17	13	19	16	22	21.1	24.5	15	20
CCMH (g/dl)	31	36	31	36	31	39	28	36	29.5	33.9	30	34
IDR (%)	14	17	17	22	14	17			12	14.5	14	17
Plaquettes ($10^3/\mu\text{l}$)	200	460	180	430	100	400	175	600	250	610	300	900
VMP (μm^3)	6.7	11.1	6.5	15	6.7	11.1			5	7	6.7	11.1