










scil Vet abc Plus ⁺

Valeurs de référence



																		
	de	à	de	à	de	à	de	à	de	à	de	à	de	à	de	à	de	à
Hémoglobine (g/dl)	15	20	8	17	9	18	11	16.5	11.5	15.1	11	16.5	11.1	18	11.1	18	8	15
Hématocrite (%)	44	57	27	47	30	50	35	50	36	47	35	50	36	52	36	52	24	46
Plaquettes (10 ³ /mm ³)	200	460	180	430	100	400	150	390	250	610	150	390	140	600	140	600	175	600
Globules rouges (10 ⁶ /mm ³)	5.5	8.5	5	10	5.5	11	3.8	5.8	5.2	6.8	3.8	5.8	5	12	5	12	5	10
Globules blancs (10 ³ /mm ³)	6	13	5	19	5	10	3.5	10	6.3	10	3.5	10	3	15	3	15	4	12
Monocytes (10 ³ /mm ³)	0	0.5	0	0.5	0	0.4												
Granulocytes (10 ³ /mm ³)	3	10	3	12	3	8												
Lymphocytes (10 ³ /mm ³)	1	3.6	1	7	1.5	4												
VGM (fl)	60	77	40	55	37	55	80	97	65	76	80	97	44	69	44	69	40	60
TCMH (pg)	17	26	13	17	13	19	26.5	33.5	21.1	24.5	26.5	33.5	12	24.5	12	24.5	16	22
CCMH (g/dl)	31	38	31	36	31	39	31.5	35	29.5	33.9	31.5	35	21.6	42	21.6	42	28	36
IDR (%)	14	17	17	22	14	20			12	15			21	25	21	25		
VMP (fl)	6.7	11.1	6.5	15	6.7	11.1			5	7			5	7	5	7		



scil animal care company
Sarl
Parc Activeum
3, rue Jacqueline Auriol
67 120 ALTORF
Tel : 03 90 20 16 40
Fax : 03 90 20 16 45
info-fr@scilvet.com
www.scilvet.com