




The scil Reflovet® Plus Chemistry Analyzer – When one or two tests are all you need.

The Reflovet® Plus Chemistry Analyzer is used in a variety of practices that need an easy to operate analyzer capable of performing a few tests quickly and accurately. Reflovet® Plus provides results for 16 parameters which allows greater flexibility for targeted test choices.

- **Small sample size - 32ul**
- **Whole blood, serum or plasma**
- **Just 3 short steps to a result**
- **Results in about 2 minutes**
- **16 low cost parameters available**
- **2 year warranty**
- **24/7 technical support**



ANALYTES



Alkaline Phosphatase, Amylase, Bilirubin, BUN, Cholesterol, Creatinine, Creatine Kinase, Glucose, GGT, GOT (AST), GPT (ALT), Hemoglobin, Potassium, Triglycerides, Uric Acid, Urea

The scil Reflovet[®] Plus Chemistry Analyzer –

Chemistry tests have never been easier, more accurate or less expensive than they are now.

JUST THREE STEPS TO A RESULT:




1

Add whole blood, plasma or serum



2

Insert the test strip



3

Results print within 2-3 minutes

TECHNICAL DATA

Dimensions:

30 cm (D) x 35 cm (L) x 21 cm (H)

Netweight:

11.68 lbs/5.3 kg

Power requirements:

115-230 V

Computer interface:

RS232 serial

Internal storage:

Up to 60 test results

Sample Volume:

32ul whole blood, serum or plasma



Quality Performance, Proven Results.