

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0724A



2024-09-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0724A-1	LOT	ME0724A-2	LOT	ME0724A-3
element HT5+	Ret# x10 ³ /μL	65.7	± 64.4	226.1	± 99.0	193.5	± 99.0
QC Mode	Ret%	1.12	± 1.12	4.82	± 2.00	8.34	± 3.50
	IRF %	10.8	± 15.0	12.3	± 13.0	10.2	± 25.0
	LFR %	89.2	± 15.0	87.7	± 15.0	89.8	± 25.0
	MFR %	8.6	± 12.0	10.9	± 10.0	8.1	± 22.0
	HFR %	2.2	± 10.0	1.4	± 10.0	2.1	± 10.0
	RHE pg	27.0	± 12.0	24.9	± 12.0	21.4	± 12.0
	IPF %	7.7	± 15.0	7.5	± 15.0	7.5	± 15.0

Before using, refer to the instruction sheet for mixing directions.