

Samsung PT10V

Profile Overview



Parameter	Electrolyte Test	Kidney Test	Large Animal Test	Liver Test	Pre-Surgical Test	Kidney Monitoring Test	Liver Monitoring Test	Diabetes Test	cCRP Test	Bovine Test	Equine Test	Comprehensive Test 17 V
A/G ratio*			X	X	X						X	X
ALB		X	X	X	X						X	X
ALP			X	X	X						X	X
ALT				X	X		X					X
AMY		X										X
AST			X	X						X	X	
B/C ratio*		X	X		X						X	X
BUN		X	X	X	X	X					X	X
CA		X								X	X	X
cCRP									X			
CHOL								X				X
CK			X							X	X	
CL	X											
CREA		X	X		X	X					X	X
FRU								X				
GGT			X	X			X				X	X
GLDH											X	
GLOB*			X	X	X						X	X
GLU		X	X	X	X			X			X	X
K	X											
LIPA												X
MG										X		
NA	X											
NA/K ratio*	X											
PHOS		X				X				X	X	X
TBIL			X	X			X				X	X
TP			X	X	X						X	X
TRIG								X				