

THE VALUES BELOW ARE FOR USE WITH THE ABAXIS PICCOLO / PICCOLO XPRESS SYSTEM

| METHOD | Abbreviation | Unit | 855UECM Level 2 | | |
|----------------------------|------------------|--------|-----------------|-----------|------------|
| | | | TARGET | LOW LIMIT | HIGH LIMIT |
| Albumin | ALB | g/dL | 4.6 | 4.2 | 5.1 |
| | | g/L | 46 | 42 | 51 |
| Alkaline phosphatase | ALP | U/L | 244 | 183 | 305 |
| Alanine aminotransferase | ALT | U/L | 96 | 74 | 118 |
| Amylase | AMY | U/L | 213 | 173 | 253 |
| Aspartate aminotransferase | AST | U/L | 92 | 71 | 113 |
| Blood urea nitrogen | BUN | mg/dL | 44 | 40 | 49 |
| | | mmol/L | 15.7 | 14.3 | 17.5 |
| Calcium | Ca | mg/dL | 11.0 | 9.8 | 12.2 |
| | | mmol/L | 2.74 | 2.45 | 3.04 |
| Cholesterol | CHOL | mg/dL | 179 | 154 | 205 |
| | | mmol/L | 4.64 | 3.99 | 5.31 |
| Creatine Kinase | CK | U/L | 528 | 422 | 634 |
| Chloride | CL ⁻ | mmol/L | 122 | 111 | 133 |
| Creatinine | CRE | mg/dL | 3.6 | 2.8 | 4.4 |
| | | µmol/L | 319 | 248 | 389 |
| Direct Bilirubin | DBIL | mg/dL | 3.2 | 2.3 | 4.0 |
| | | µmol/L | 54.7 | 39.3 | 68.4 |
| Gamma GT | GGT | U/L | 130 | 102 | 159 |
| Glucose | GLU | mg/dL | 248 | 208 | 288 |
| | | mmol/L | 13.8 | 11.5 | 16.0 |
| High-Density Lipoprotein | HDL | mg/dL | 69 | 59 | 79 |
| | | mmol/L | 1.8 | 1.5 | 2.1 |
| Potassium | K ⁺ | mmol/L | 6.1 | 5.5 | 6.7 |
| Lactate Dehydrogenase | LDH | U/L | 664 | 545 | 784 |
| Magnesium | Mg | mg/dL | 4.1 | 3.5 | 4.8 |
| | | mmol/L | 1.7 | 1.4 | 2.0 |
| Sodium | Na ⁺ | mmol/L | 163 | 153 | 173 |
| Phosphorus | PHOS | mg/dL | 7.0 | 5.7 | 8.2 |
| | | mmol/L | 2.3 | 1.8 | 2.7 |
| Total Bilirubin | TBIL | mg/dL | 4.9 | 3.6 | 6.2 |
| | | µmol/L | 83.8 | 61.5 | 106 |
| Bicarbonate | TCO ₂ | mmol/L | 24 | 18 | 30 |
| Total Protein | TP | g/dL | 6.4 | 5.7 | 7.1 |
| | | g/L | 64 | 57 | 71 |
| Triglycerides | TRIG | mg/dL | 273 | 224 | 322 |
| | | mmol/L | 3.1 | 2.5 | 3.6 |
| Uric Acid | UA | mg/dL | 9.6 | 8.3 | 10.9 |
| | | µmol/L | 571 | 494 | 649 |
| Lactate | LAC | mmol/L | 4.2 | 3.4 | 4.9 |
| | | mg/dL | 37.8 | 30.6 | 44.1 |
| C-Reactive Protein | CRP | mg/L | 25.5 | 18.6 | 32.4 |
| | | mg/dL | 2.6 | 1.9 | 3.2 |

13 Jul '16 ne