

HORIBA AMX 335 L
HORIBA ABX SAS Rev 1

CONTROL L

(Exp.)

2022-08-05
(YYYY-MM-DD)

IVD

		LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
	UNITES UNITS	Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.7	1.7						± 0.4	
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.26	2.24						± 0.15	
HB HGB	g/dl	5.8	5.7						± 0.4	
	g/l	58	57						± 4	
	mmol/l	3.60	3.54						± 0.25	
HT HCT	%	16.0	15.9						± 2.0	
	l/l	0.160	0.159						± 0.020	
VGM MCV	µm ³ ; fl	71	71						± 4	
TGMH MCH	pg	25.7	25.4						± 2.0	
	fmol	1.59	1.58						± 0.12	
CCMH MCHC	g/dl	36.1	35.8						± 3.0	
	g/l	361	358						± 30	
	mmol/l	22.45	22.26						± 1.86	
IDR RDW	%	N/A	15.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	71	75						± 20	
VPM MPV	µm ³ ; fl	N/A	9.1						± 2.0	
LYMPHO	#	N/A	0.93						± 0.40	
	%	N/A	54.5						± 8.0	
MONO	#	N/A	0.21						± 0.20	
	%	N/A	12.5						± 6.0	
GRANULO	#	N/A	0.56						± 0.40	
	%	N/A	33.0						± 7.0	

LOT MX 435 N
Rev 1

CONTROL N

(Exp.)

2022-08-05
(YYYY-MM-DD)

		LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
PARAMETRES PARAMETERS	UNITES UNITS	Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.3						± 0.8	
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.63	4.63						± 0.18	
HB HGB	g/dl	13.4	13.4						± 0.5	
	g/l	134	134						± 5	
	mmol/l	8.32	8.32						± 0.31	
HT HCT	%	37.0	37.5						± 2.5	
	l/l	0.370	0.375						± 0.025	
VGM MCV	µm ³ ; fl	80	81						± 4	
TGMH MCH	pg	28.9	28.9						± 2.0	
	fmol	1.80	1.80						± 0.12	
CCMH MCHC	g/dl	36.2	35.7						± 3.0	
	g/l	362	357						± 30	
	mmol/l	22.47	22.19						± 1.86	



HORIBA ABX SAS

HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE



TEMP	33.07	1.07	14.15	16															± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		245	255															± 40
VPM MPV	μm ³ ; fl		N/A	8.8															± 2.0
LYMPHO	#		N/A	2.04															± 0.40
	%		N/A	28.0															± 5.0
MONO	#		N/A	0.62															± 0.30
	%		N/A	8.5															± 4.0
GRANULO	#		N/A	4.64															± 0.60
	%		N/A	63.5															± 6.0

Ref: TEMP-0830 Rev:41 BACK

LOT	MX 435 H Rev 1	CONTROL	H	(Exp.) 2022-08-05 (YYYY-MM-DD)
------------	-------------------	----------------	---	--

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet												TOLERANCES TOLERANCE
		MICROS						MICROS						
		Sci Vet ABC			Sci Vet ABC Plus ESV 60									
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	18.8			19.0									± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.65			5.65									± 0.20
HB HGB	g/dl	18.1			18.0									± 0.6
	g/l	181			180									± 6
	mmol/l	11.24			11.18									± 0.37
HT HCT	%	50.3			50.3									± 3.0
	l/l	0.503			0.503									± 0.030
VGM MCV	μm ³ ; fl	89			89									± 4
TGMH MCH	pg	32.0			31.9									± 2.0
	fmol	1.99			1.98									± 0.12
CCMH MCHC	g/dl	36.0			35.8									± 3.0
	g/l	360			358									± 30
	mmol/l	22.35			22.23									± 1.86
IDR RDW	%	N/A			13.5									± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	500			525									± 45
VPM MPV	μm ³ ; fl	N/A			8.5									± 2.0
LYMPHO	#	N/A			2.09									± 1.00
	%	N/A			11.0									± 5.0
MONO	#	N/A			0.95									± 0.80
	%	N/A			5.0									± 4.0
GRANULO	#	N/A			15.96									± 1.40
	%	N/A			84.0									± 7.0