

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MX 447 L
Rev 1

CONTROL L

 (Exp.)

2024-08-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet							TOLERANCES TOLERANCE
		MICROS							
		Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27							± 0.15
HB HGB	g/dl	5.9							± 0.4
	g/l	59							± 4
	mmol/l	3.66							± 0.25
HT HCT	%	16.3							± 2.0
	l/l	0.163							± 0.020
VGM MCV	µm ³ ; fl	72							± 4
TGMH MCH	pg	26.0							± 2.0
	fmol	1.61							± 0.12
CCMH MCHC	g/dl	36.1							± 3.0
	g/l	361							± 30
	mmol/l	22.42							± 1.86
IDR RDW *	%	15.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	74							± 20
VPM MPV *	µm ³ ; fl	8.0							± 2.0
LYMPHO *	#	0.96							± 0.40
	%	53.5							± 8.0
MONO *	#	0.23							± 0.20
	%	13.0							± 6.0
GRANULO *	#	0.60							± 0.40
	%	33.5							± 7.0

* N/A for Micros Scil Vet ABC

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17



LOT

MX 447 N
Rev 1

CONTROL

N

(Exp.)

2024-08-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.1							± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l		4.54							± 0.18
HB HGB	g/dl		13.2							± 0.5
	g/l		132							± 5
	mmol/l		8.20							± 0.31
HT HCT	%		36.8							± 2.5
	l/l		0.368							± 0.025
VGM MCV	µm ³ ; fl		81							± 4
TGMH MCH	pg		29.1							± 2.0
	fmol		1.81							± 0.12
CCMH MCHC	g/dl		35.9							± 3.0
	g/l		359							± 30
	mmol/l		22.29							± 1.86
IDR RDW *	%		14.5							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		262							± 40
VPM MPV *	µm ³ ; fl		7.5							± 2.0
LYMPHO *	#		2.02							± 0.40
	%		28.5							± 5.0
MONO *	#		0.71							± 0.45
	%		10.0							± 6.0
GRANULO *	#		4.37							± 0.60
	%		61.5							± 6.0

* N/A for Micros Scil Vet ABC

1300104143
BACK / VERSO
Ref: TEMP-0630 Rev.43

LOT MX 447 H
Rev 1

CONTROL H

(Exp.)

2024-08-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		19.7							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l		5.73							± 0.20
HB HGB	g/dl		18.1							± 0.6
	g/l		181							± 6
	mmol/l		11.24							± 0.37
HT HCT	%		50.4							± 3.0
	l/l		0.504							± 0.030
VGM MCV	µm ³ ; fl		88							± 4
TGMH MCH	pg		31.6							± 2.0
	fmol		1.96							± 0.12
CCMH MCHC	g/dl		35.9							± 3.0
	g/l		359							± 30
	mmol/l		22.29							± 1.86
IDR RDW *	%		13.5							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		498							± 45
VPM MPV *	µm ³ ; fl		7.4							± 2.0
LYMPHO *	#		2.17							± 1.00
	%		11.0							± 5.0
MONO *	#		1.38							± 0.80
	%		7.0							± 4.0
GRANULO *	#		16.15							± 1.40
	%		82.0							± 7.0

* N/A for Micros Scil Vet ABC