

# Abaxis VetScan®

## Profiles and Parameters



Parameter	Mammalian Liver	Equine Plus	Critical Care Plus	Large Animal	Comprehensive Diagnostic	Prep Profile II	Thyroxine / Cholesterol	Avian / Reptilian Plus	Kidney Profile Plus	Electrolyte Plus	Preventive Care Profile Plus	Phenobarbital
ALB	x	x		x	x			x	x		x	x
ALP	x			x	x	x					x	x
AMYL					x							
BUN	x	x	x	x	x	x			x		x	x
CA		x		x	x			x	x		x	
CHOL	x						x					
CK		x		x				x				
Cl			x						x	x	x	
CREA		x	x		x	x			x		x	
Bile Acid	x							x				
GGT	x	x		x								x
GLOB		x		x	x			x			x	
GLU		x	x		x	x		x	x		x	
GOT / AST		x		x				x			x	x
GPT / ALT	x		x		x	x					x	x
K		x	x		x			x	x	x	x	
MG				x								
NA		x	x		x			x	x	x	x	
PHB												x
PHOS				x	x			x	x			
T4							x					
TBIL	x	x			x						x	x
tCO <sub>2</sub>		x	x						x	x	x	
TP		x		x	x	x		x			x	
UA								x				