

Radox Quality Control Human Level 2 (5ml)

For the control for clinical chemistry parameters of Element RC

FOR VETERINARY USE ONLY!

Antech[™]



Itemcode 111843

LOT 1750UN

Expiry date 28.06.2028

Handling

To be reconstituted with 5 ml of distilled water, mix the control gently and wait for 30min.

Analyte	Unit	Target	MIN	MAX	Unit	Target	MIN	MAX
ALB	g/l	37,6	33,8	41,4	g/dl	3,76	3,38	4,14
ALP	U/L	178	160	196	U/L	178	160	196
ALT	U/L	54	46	62	U/L	54	46	62
AMY	U/L	76	68	84	U/L	76	68	84
AST	U/L	24	20	28	U/L	24	20	28
BUN	mmol/l	6,32	5,69	6,95	mg/dl	17,7	15,93	19,46
Ca	mmol/l	2,15	1,94	2,37	mg/dl	8,62	7,78	9,50
CK	U/L	210	189	231	U/L	210	189	231
Cl	mmol/l	96,3	86,7	105,9	mmol/l	96,3	86,7	105,9
Crea	μmol/l	103,5	88,0	119,0	mg/dl	1,18	1,0	1,35
GGT	U/L	47	42	52	U/L	47	42	52
GLU	mmol/l	5,94	5,35	6,53	mg/dl	107	96	118
K	mmol/l	4,42	3,76	5,08	mmol/l	4,42	3,76	5,08
LAC	mmol/l	1,71	1,45	1,97	mmol/l	1,71	1,45	1,97
LDH	U/L	232	209	255	U/L	232	209	255
Mg	mmol/l	0,95	0,86	1,05	mg/dl	2,31	2,09	2,55
Na	mmol/l	146,3	139	153,6	mmol/l	146,3	139	153,6
PHOS	mmol/l	1,56	1,4	1,72	mg/dl	4,83	4,34	5,33
TB	μmol/l	30,7	27,6	33,8	mg/dl	1,80	1,61	1,98
TBA	μmol/l	28,63	24,34	32,92	μmol/l	28,63	24,34	32,92
TC	mmol/l	3,74	3,37	4,11	mg/dl	145	130	159
tCO2	mmol/l	14	11,9	16,1	mmol/l	14	11,9	16,1
TG	mmol/l	1,23	1,11	1,35	mg/dl	108,85	98,23	119,47
TP	g/l	52,3	47,1	57,5	g/dl	5,23	4,71	5,75
UA	μmol/l	326,48	277,51	375,45	μmol/l	326,48	277,51	375,45