

scil vCell 5 QC

Control materials for scil vCell 5
veterinary hematology analyzers

CONTROL

LOT N5P129



05-03-2026

Software version: 1.6.2017.0

Assay Sheet revision: 05-12-2025

Parameters	Units	LOW				NORMAL				HIGH			
		mean	limit	min	max	mean	limit	min	max	mean	limit	min	max
WBC / GB	10 ³ /μL & 10 ⁹ /L	3.5	± 0.6	2.9	4.1	7.3	± 0.8	6.5	8.1	20.1	± 2.2	17.9	22.3
LYM%	%	42.6	± 7.0	35.6	49.6	20.3	± 6.0	14.3	26.3	11.4	± 6.0	5.4	17.4
MON%	%	7.0	± 5.0	2.0	12.0	4.8	± 4.5	0.3	9.3	4.2	± 3.0	1.2	7.2
NEU%	%	46.0	± 7.0	39.0	53.0	69.1	± 7.0	62.1	76.1	78.1	± 8.0	70.1	86.1
EOS%	%	3.1	± 2.5	0.6	5.6	5.0	± 4.0	1.0	9.0	5.7	± 4.5	1.2	10.2
BAS%	%	1.3	± 0.7	0.6	2.0	0.8	± 0.5	0.3	1.3	0.6	± 0.3	0.3	0.9
LYM#	10 ³ /μL & 10 ⁹ /L	1.5	± 0.4	1.1	1.9	1.5	± 0.6	0.9	2.1	2.3	± 1.3	1.0	3.6
MON#	10 ³ /μL & 10 ⁹ /L	0.2	± 0.2	0.0	0.4	0.4	± 0.4	0.0	0.8	0.8	± 0.8	0.0	1.6
NEU#	10 ³ /μL & 10 ⁹ /L	1.6	± 0.5	1.1	2.1	4.9	± 1.3	3.6	6.2	15.8	± 3.0	12.8	18.8
EOS#	10 ³ /μL & 10 ⁹ /L	0.1	± 0.1	0.0	0.2	0.4	± 0.3	0.1	0.7	1.1	± 1.0	0.1	2.1
BAS#	10 ³ /μL & 10 ⁹ /L	0.1	± 0.1	0.0	0.2	0.1	± 0.1	0.0	0.2	0.1	± 0.1	0.0	0.2
RBC / GR	10 ⁶ /μL & 10 ¹² /L	2.06	± 0.20	1.86	2.26	4.53	± 0.32	4.21	4.85	5.06	± 0.40	4.66	5.46
HGB	g/dL	5.5	± 0.4	5.1	5.9	12.7	± 0.6	12.1	13.3	15.5	± 0.8	14.7	16.3
	g/L	55	± 4	51	59	127	± 6	121	133	155	± 8	147	163
	mmol/L	3.41	± 0.25	3.16	3.66	7.88	± 0.37	7.51	8.25	9.62	± 0.50	9.12	10.12
HCT	%	18.5	± 2.7	15.8	21.2	45.3	± 5.4	39.9	50.7	52.6	± 6.3	46.3	58.9
	L/L	0.19	± 0.03	0.16	0.22	0.45	± 0.06	0.39	0.51	0.53	± 0.07	0.46	0.60
MCV / VGM	fL	90	± 5	85	95	100	± 5	95	105	104	± 5	99	109
MCH / TCMH	pg	26.7	± 3.8	22.9	30.5	28.0	± 2.8	25.2	30.8	30.6	± 3.0	27.6	33.6
	fmol	1.66	± 0.24	1.42	1.90	1.74	± 0.17	1.57	1.91	1.90	± 0.19	1.71	2.09
MCHC / CCMH	g/dL	29.7	± 4.5	25.2	34.2	28.0	± 3.7	24.3	31.7	29.5	± 3.6	25.9	33.1
	g/L	297	± 45	252	342	280	± 37	243	317	295	± 36	259	331
	mmol/L	18.4	± 2.8	15.6	21.2	17.4	± 2.3	15.1	19.7	18.3	± 2.2	16.1	20.5
RDWcv / IDR	%	17.0	± 3.0	14.0	20.0	15.5	± 2.8	12.7	18.3	15.2	± 2.5	12.7	17.7
PLT	10 ³ /μL & 10 ⁹ /L	62	± 25	37	87	214	± 50	164	264	472	± 70	402	542
PCT / Tct	%	0.05	± 0.04	0.01	0.09	0.17	± 0.08	0.09	0.25	0.38	± 0.15	0.23	0.53
MPV / VPM	fL	7.9	± 2.0	5.9	9.9	8.0	± 2.0	6.0	10.0	8.0	± 2.0	6.0	10.0
PDWcv / IDP	%	57.0	± 9.0	48.0	66.0	58.0	± 7.0	51.0	65.0	58.0	± 6.0	52.0	64.0

NNEU
xy
72/142

LOW



NORMAL



HIGH



How to use the QR codes:

1. Start the analyzer, wait for the Main Menu.
2. Go to Daily Routine.
3. Tap menu (≡) button in lower right corner.
4. Tap "Read QR".
5. Align the code on the screen so that only one is visible entirely, aligned parallel with the camera and the front panel.
6. The analyzer will acknowledge successful scanning with a message.
7. Repeat the process for all three levels.

Scanning a QR code multiple times will NOT create multiple QC bank entries.

How to upload QRC files:

1. Copy QRC files to the root folder of a USB stick.
2. Connect the stick with the analyzer ON.
3. Go to Daily Routine and tap menu (≡) button in lower right corner.
4. Select "Load QR". Successful loading will be acknowledged by a message.



For further information, please refer to the instructions for use.



scil animal care company GmbH
Dina-Weissmann-Allee
668519 Viernheim
GERMANY

