

LOT		MX 455 L		CONTROL		L		 (Exp.)		2025-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		1.8								± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.23								± 0.15
HB HGB	g/dl		5.8								± 0.4
	g/l		58								± 4
	mmol/l		3.60								± 0.25
HT HCT	%		15.8								± 2.0
	l/l		0.158								± 0.020
VGM MCV	µm ³ ; fl		71								± 4
TGMH MCH	pg		26.0								± 2.0
	fmol		1.62								± 0.12
CCMH MCHC	g/dl		36.6								± 3.0
	g/l		366								± 30
	mmol/l		22.75								± 1.86
IDR RDW *	%		15.5								± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		74								± 20
VPM MPV *	µm ³ ; fl		8.3								± 2.0
LYMPHO *	#		0.98								± 0.40
	%		54.5								± 8.0
MONO *	#		0.23								± 0.21
	%		13.0								± 11.7
GRANULO *	#		0.59								± 0.40
	%		32.5								± 10.0

* N/A for Micros Scil Vet ABC

LOT		MX 455 N		CONTROL		N		(Exp.)		2025-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.6								± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.55								± 0.18
HB HGB	g/dl		13.4								± 0.5
	g/l		134								± 5
	mmol/l		8.32								± 0.31
HT HCT	%		37.3								± 2.5
	l/l		0.373								± 0.025
VGM MCV	µm ³ ; fl		82								± 4
TGMH MCH	pg		29.5								± 2.0
	fmol		1.83								± 0.12
CCMH MCHC	g/dl		35.9								± 3.0
	g/l		359								± 30
	mmol/l		22.30								± 1.86
IDR RDW *	%		15.0								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		266								± 40
VPM MPV *	µm ³ ; fl		7.8								± 2.0
LYMPHO *	#		2.01								± 0.40
	%		26.5								± 5.0
MONO *	#		0.76								± 0.68
	%		10.0								± 9.0
GRANULO *	#		4.83								± 0.76
	%		63.5								± 10.0

* N/A for Micros Scil Vet ABC

1300104143
BACK / VERSO
Ref: TEMP-0830 Rev.43

LOT		MX 455 H		CONTROL		H		 (Exp.)		2025-12-05 (YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		19.8								± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.69								± 0.20
HB HGB	g/dl		17.9								± 0.6
	g/l		179								± 6
	mmol/l		11.12								± 0.37
HT HCT	%		50.6								± 3.0
	l/l		0.506								± 0.030
VGM MCV	µm ³ ; fl		89								± 4
TGMH MCH	pg		31.5								± 2.0
	fmol		1.95								± 0.12
CCMH MCHC	g/dl		35.3								± 3.0
	g/l		353								± 30
	mmol/l		21.95								± 1.86
IDR RDW *	%		14.0								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		512								± 45
VPM MPV *	µm ³ ; fl		7.6								± 2.0
LYMPHO *	#		2.28								± 1.00
	%		11.5								± 5.0
MONO *	#		1.19								± 1.07
	%		6.0								± 5.4
GRANULO *	#		16.34								± 1.98
	%		82.5								± 10.0

* N/A for Micros Scil Vet ABC