

LOT		MX 456 L		CONTROL		L		 (Exp.)		2026-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		1.9								± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.28								± 0.15
HB HGB	g/dl		6.0								± 0.4
	g/l		60								± 4
	mmol/l		3.73								± 0.25
HT HCT	%		16.4								± 2.0
	l/l		0.164								± 0.020
VGM MCV	µm ³ ; fl		72								± 4
TGMH MCH	pg		26.3								± 2.0
	fmol		1.63								± 0.12
CCMH MCHC	g/dl		36.5								± 3.0
	g/l		365								± 30
	mmol/l		22.70								± 1.86
IDR RDW *	%		15.0								± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		75								± 20
VPM MPV *	µm ³ ; fl		7.2								± 2.0
LYMPHO *	#		0.99								± 0.40
	%		52.0								± 8.0
MONO *	#		0.26								± 0.23
	%		13.5								± 12.2
GRANULO *	#		0.66								± 0.40
	%		34.5								± 10.0

* N/A for Micros Scil Vet ABC

LOT		MX 456 N		CONTROL		N		(Exp.)		2026-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.2								± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.56								± 0.18
HB HGB	g/dl		13.4								± 0.5
	g/l		134								± 5
	mmol/l		8.32								± 0.31
HT HCT	%		37.4								± 2.5
	l/l		0.374								± 0.025
VGM MCV	µm ³ ; fl		82								± 4
TGMH MCH	pg		29.4								± 2.0
	fmol		1.82								± 0.12
CCMH MCHC	g/dl		35.8								± 3.0
	g/l		358								± 30
	mmol/l		22.25								± 1.86
IDR RDW *	%		14.0								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		244								± 40
VPM MPV *	µm ³ ; fl		7.3								± 2.0
LYMPHO *	#		2.05								± 0.40
	%		28.5								± 5.0
MONO *	#		0.58								± 0.52
	%		8.0								± 7.2
GRANULO *	#		4.57								± 0.72
	%		63.5								± 10.0

* N/A for Micros Scil Vet ABC

1300.104143
BACK / VERSO
Ref: TEMP-0830 Rev.43

LOT		MX 456 H		CONTROL		H		 (Exp.)		2026-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.4								± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.72								± 0.20
HB HGB	g/dl		18.3								± 0.6
	g/l		183								± 6
	mmol/l		11.36								± 0.37
HT HCT	%		50.9								± 3.0
	l/l		0.509								± 0.030
VGM MCV	µm ³ ; fl		89								± 4
TGMH MCH	pg		32.0								± 2.0
	fmol		1.99								± 0.12
CCMH MCHC	g/dl		35.9								± 3.0
	g/l		359								± 30
	mmol/l		22.32								± 1.86
IDR RDW *	%		13.5								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		526								± 45
VPM MPV *	µm ³ ; fl		7.2								± 2.0
LYMPHO *	#		2.24								± 1.00
	%		11.0								± 5.0
MONO *	#		0.92								± 0.83
	%		4.5								± 4.1
GRANULO *	#		17.24								± 2.04
	%		84.5								± 10.0

* N/A for Micros Scil Vet ABC