

LOT		MX 457 L		CONTROL		L		 (Exp.)		2026-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		2.0								± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.24								± 0.15
HB HGB	g/dl		5.8								± 0.4
	g/l		58								± 4
	mmol/l		3.60								± 0.25
HT HCT	%		15.9								± 2.0
	l/l		0.159								± 0.020
VGM MCV	µm ³ ; fl		71								± 4
TGMH MCH	pg		25.9								± 2.0
	fmol		1.61								± 0.12
CCMH MCHC	g/dl		36.5								± 3.0
	g/l		365								± 30
	mmol/l		22.65								± 1.86
IDR RDW *	%		15.0								± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		76								± 20
VPM MPV *	µm ³ ; fl		8.3								± 2.0
LYMPHO *	#		1.11								± 0.40
	%		55.5								± 8.0
MONO *	#		0.26								± 0.23
	%		13.0								± 11.7
GRANULO *	#		0.63								± 0.40
	%		31.5								± 10.0

* N/A for Micros Scil Vet ABC

LOT		MX 457 N		CONTROL		N		(Exp.)		2026-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.4								± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.60								± 0.18
HB HGB	g/dl		13.4								± 0.5
	g/l		134								± 5
	mmol/l		8.32								± 0.31
HT HCT	%		37.3								± 2.5
	l/l		0.373								± 0.025
VGM MCV	µm ³ ; fl		81								± 4
TGMH MCH	pg		29.1								± 2.0
	fmol		1.81								± 0.12
CCMH MCHC	g/dl		36.0								± 3.0
	g/l		360								± 30
	mmol/l		22.33								± 1.86
IDR RDW *	%		15.0								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		280								± 40
VPM MPV *	µm ³ ; fl		8.1								± 2.0
LYMPHO *	#		1.89								± 0.40
	%		25.5								± 5.0
MONO *	#		0.74								± 0.67
	%		10.0								± 9.0
GRANULO *	#		4.77								± 0.74
	%		64.5								± 10.0

* N/A for Micros Scil Vet ABC

1300104143
BACK / VERSO
Ref: TEMP-0830 Rev.43

LOT		MX 457 H		CONTROL		H		 (Exp.)		2026-04-05 (YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
			MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.3								± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.69								± 0.20
HB HGB	g/dl		18.0								± 0.6
	g/l		180								± 6
	mmol/l		11.18								± 0.37
HT HCT	%		50.1								± 3.0
	l/l		0.501								± 0.030
VGM MCV	µm ³ ; fl		88								± 4
TGMH MCH	pg		31.6								± 2.0
	fmol		1.96								± 0.12
CCMH MCHC	g/dl		35.9								± 3.0
	g/l		359								± 30
	mmol/l		22.32								± 1.86
IDR RDW *	%		14.0								± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		527								± 45
VPM MPV *	µm ³ ; fl		7.8								± 2.0
LYMPHO *	#		2.13								± 1.00
	%		10.5								± 5.0
MONO *	#		1.12								± 1.01
	%		5.5								± 5.0
GRANULO *	#		17.05								± 2.03
	%		84.0								± 10.0

* N/A for Micros Scil Vet ABC